

STERILE CONTAINER SYSTEMS





Standard Model

» 1/1 Full Size

1/1 STANDARD MODEL

W: 285 mm

L: 600 mm




Specifications		Lid (Perforated)	Lid (Perforated)	Safety Lid	Safety Lid
Height					
		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)	Box (Perforated)
H:110 mm		Z6-010	Z6-040	Z6-150	Z6-180
		Z6-011	Z6-041	Z6-151	Z6-181
		Z6-012	Z6-042	Z6-152	Z6-182
		Z6-013	Z6-043	Z6-153	Z6-183
		Z6-014	Z6-044	Z6-154	Z6-184
		Z6-015	Z6-045	Z6-155	Z6-185
H:135 mm		Z6-016	Z6-046	Z6-156	Z6-186
		Z6-017	Z6-047	Z6-157	Z6-187
		Z6-018	Z6-048	Z6-158	Z6-188
		Z6-019	Z6-049	Z6-159	Z6-189
		Z6-020	Z6-050	Z6-160	Z6-190
		Z6-021	Z6-051	Z6-161	Z6-191
H:155 mm		Z6-022	Z6-052	Z6-162	Z6-192
		Z6-023	Z6-053	Z6-163	Z6-193
		Z6-024	Z6-054	Z6-164	Z6-194
		Z6-025	Z6-055	Z6-165	Z6-195
		Z6-026	Z6-056	Z6-166	Z6-196
		Z6-027	Z6-057	Z6-167	Z6-197
H:205 mm		Z6-028	Z6-058	Z6-168	Z6-198
		Z6-029	Z6-059	Z6-169	Z6-199
		Z6-030	Z6-060	Z6-170	Z6-200
		Z6-031	Z6-061	Z6-171	Z6-201
		Z6-032	Z6-062	Z6-172	Z6-202
		Z6-033	Z6-063	Z6-173	Z6-203
H:255 mm		Z6-034	Z6-064	Z6-174	Z6-204
		Z6-035	Z6-065	Z6-175	Z6-205
		Z6-036	Z6-066	Z6-176	Z6-206
		Z6-037	Z6-067	Z6-177	Z6-207
		Z6-038	Z6-068	Z6-178	Z6-208
		Z6-039	Z6-069	Z6-179	Z6-209


Ps: All container measurements cover the outer size of the lids.

BOX (perforated)


outside dimension



outside dimension




inside dimension



code	outside dimension	inside dimension	weight
Z6-100	600 x 272 h:100 mm	543 x 258 h:83 mm	2216 g
Z6-101	600 x 272 h:116 mm	543 x 258 h:99 mm	2362 g
Z6-102	600 x 272 h:138 mm	543 x 258 h:121 mm	2520 g
Z6-103	600 x 272 h:190 mm	543 x 258 h:173 mm	2873 g
Z6-104	600 x 272 h:235 mm	543 x 258 h:218 mm	3290 g

BOX (non perforated)

outside dimension



outside dimension



inside dimension



code	outside dimension	inside dimension	weight
Z6-105	600 x 272 h:100 mm	543 x 258 h:94 mm	1900 g
Z6-106	600 x 272 h:116 mm	543 x 258 h:110 mm	2096 g
Z6-107	600 x 272 h:138 mm	543 x 258 h:132 mm	2205 g
Z6-108	600 x 272 h:190 mm	543 x 258 h:184 mm	2648 g
Z6-109	600 x 272 h:235 mm	543 x 258 h:229 mm	3040 g

PERFORATED LIDS

code	color	weight
Z6-120		1310 g
Z6-121		
Z6-122		
Z6-123		
Z6-124		
Z6-125		

Dimension: 602 x 285 mm

REMOVAL LIDS (non perforated)

code	color	weight
Z6-130		1046 g
Z6-131		
Z6-132		
Z6-133		
Z6-134		
Z6-135		

Dimension: 602 x 285 mm


SAFETY LIDS

code	color	weight
Z6-230		988 g
Z6-231		
Z6-232		
Z6-233		
Z6-234		
Z6-235		

Dimension: 606 x 288 mm

BASKET MADE FROM STAINLESS STEEL WIRE MESH


WIRE BASKET



code	outside dimension
Z2-460	540 x 245 h:30 mm
Z2-461	540 x 245 h:50 mm
Z2-462	540 x 245 h:70 mm
Z2-463	540 x 245 h:100 mm
Z2-450	480 x 245 h:30 mm
Z2-451	480 x 245 h:50 mm
Z2-452	480 x 245 h:70 mm
Z2-453	480 x 245 h:100 mm



WIRE BASKET WITH LID




code	outside dimension
Z3-1010	540 x 245 h:35 mm
Z3-1011	540 x 245 h:55 mm
Z3-1012	540 x 245 h:75 mm
Z3-1013	540 x 245 h:105 mm
Z3-1014	ONLY LID, for baskets
Z3-1020	480 x 245 h:35 mm
Z3-1021	480 x 245 h:55 mm
Z3-1022	480 x 245 h:75 mm
Z3-1023	480 x 245 h:105 mm
Z3-1024	ONLY LID, for baskets




WIRE BASKET WITH PERFORATED METAL PLATE

WIRE BASKET



code	outside dimension
Z2-461 P	540 x 245 h:50 mm
Z2-462 P	540 x 245 h:70 mm
Z2-463 P	540 x 245 h:100 mm
Z2-451 P	480 x 245 h:50 mm
Z2-452 P	480 x 245 h:70 mm
Z2-453 P	480 x 245 h:100 mm



WIRE BASKET WITH LID




code	outside dimension
Z3-1011 P	540 x 245 h:55 mm
Z3-1012 P	540 x 245 h:75 mm
Z3-1013 P	540 x 245 h:105 mm
Z3-1014 P	ONLY LID, for baskets
Z3-1021 P	480 x 245 h:55 mm
Z3-1022 P	480 x 245 h:75 mm
Z3-1023 P	480 x 245 h:105 mm
Z3-1024 P	ONLY LID, for baskets





BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
Box and Lid 	Z3-720	546 x 259 h:44 mm
	Z3-721	546 x 259 h:64 mm
	Z3-722	546 x 259 h:84 mm
	Z3-723	546 x 259 h:114 mm
Box 	Z3-740	540 x 253 h:40 mm
	Z3-741	540 x 253 h:60 mm
	Z3-742	540 x 253 h:80 mm
	Z3-743	540 x 253 h:110 mm
Lid 	Z3-760	ONLY LID, for baskets



BASKET (with stand)



	code	outside dimension
Box and Lid 	Z3-724	491 x 259 h:44 mm
	Z3-725	491 x 259 h:64 mm
	Z3-726	491 x 259 h:84 mm
	Z3-727	491 x 259 h:114 mm
Box 	Z3-744	485 x 253 h:40 mm
	Z3-745	485 x 253 h:60 mm
	Z3-746	485 x 253 h:80 mm
	Z3-747	485 x 253 h:110 mm
Lid 	Z3-761	ONLY LID, for baskets

STAND



code	height
Z3-767	20 mm



STAND



code	height
Z3-768	10 mm



STAND



code	height
Z3-766	5 mm



IDENTIFICATION (aluminium)

LABEL HOLDER



code	color
Z2-560	Red
Z2-561	Blue
Z2-562	Green
Z2-563	Yellow
Z2-564	Grey
Z2-565	Black

LABELS



code	color
Z2-570	Red
Z2-571	Blue
Z2-572	Green
Z2-573	Yellow
Z2-574	Grey
Z2-575	Black

Coloured identification labels
Engraved to your specification



- ➡ Grey coloured recommended for use in Sterrad.
- Color fading By effect of Hydrogen Peroxide Gas Plasma on other colors.

DIVIDERS



Z2-428
540 x 250 x 130 mm
Z2-429
540 x 250 x 180 mm

WIRE BASE


Z2-440
545 x 255 x 30 mm

DIVIDERS


Z2-425
540 x 250 x 130 mm
Z2-426
540 x 250 x 180 mm

RACK (for endoscopy)


Z3-1090
505 x 250 x 125 mm
• Holds 10 endoscopic instruments
• Fits into a 1/1 size 20 cm deep container
• Wire basket 1/1 50 mm
• Container 1/1 200 mm
• Silicone mat 520 x 230 mm

DRAPE RETAINER


Z2-420
545 x 255 x 30 mm

DRAPE RETAINER


Z2-400
545 x 255 mm

ALUMINIUM BASKET (with lid)



Z3-1119
245 x 130 x 55 mm

ALUMINIUM BASKET (with lid)



Z3-1120
245 x 130 x 70 mm

ALUMINIUM BASKET (with lid)


Z3-1121
245 x 130 x 105 mm

SILICONE MAT


Z2-526
520 x 230 mm
Z2-528
475 x 230 mm

SILICONE MESH MAT


Z2-531
485 x 250 mm
Z2-533
534 x 247 mm

SILICONE HANDLE



Z2-710
Z2-711
Z2-712
Z2-713
Z2-714
Z2-715

PAPER FILTER


Z2-740
Single use 100 pcs with indicator. Insidemed mark.
• Diameter: 184 mm

TEFLON FILTER



Z2-742
Number of sterilization cycle > 2000
• Diameter: 194 mm

FILTER DISPENSER


Z2-730

IDENTIFICATION LABELS (aluminium)


Z2-700
Z2-701
Z2-702
Z2-703
Z2-704
Z2-705

LABELS (paper)


Z2-673
100 pcs. With indicator.
For small size and special containers


LABELS (paper)


Z2-672
100 pcs. With indicator.
For full and Half size containers

(-)

AUTO LOCK MECHANISM

Z1-900
(2 Pcs)

EMULATING INDICATOR CARD (class 6)


Z2-674
Paper 250 Pcs.


SECURITY SEAL (blue)

Z2-663
100 pcs. with indicator. made of plastic


SECURITY SEAL (blue)

Z2-664
100 pcs. made of plastic


SECURITY SEAL (blue)

Z2-661
100 pcs. with number. made of plastic


SECURITY SEAL (red)

Z2-662
100 pcs. with indicator. made of plastic


SECURITY SEAL (red)

Z2-660
100 pcs. made of plastic


SECURITY SEAL (red)

Z2-665
100 pcs. with number. made of plastic


DISPOSAL STAMP

Z2-717
40 x 35 x 9 mm
Requested number or letter must be informed during the order.


LABEL HOLDER

Z2-716
50 x 35 x 13 mm
Holder for identification labels

(-)



STERILE CONTAINER SYSTEMS





Standard Model




» 3/4 Quarter Size

3/4 STANDARD MODEL

W: 285 mm

L: 475 mm




Specifications		Lid (Perforated)		Safety Lid	
					
Height		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)	Box (Perforated)
H:110 mm		Z6-240 	Z6-260 	Z6-350 	Z6-370 
		Z6-241 	Z6-261 	Z6-351 	Z6-371 
		Z6-242 	Z6-262 	Z6-352 	Z6-372 
		Z6-243 	Z6-263 	Z6-353 	Z6-373 
		Z6-244 	Z6-264 	Z6-354 	Z6-374 
		Z6-245 	Z6-265 	Z6-355 	Z6-375 
H:130 mm		Z6-246 	Z6-266 	Z6-356 	Z6-376 
		Z6-247 	Z6-267 	Z6-357 	Z6-377 
		Z6-248 	Z6-268 	Z6-358 	Z6-378 
		Z6-249 	Z6-269 	Z6-359 	Z6-379 
		Z6-250 	Z6-270 	Z6-360 	Z6-380 
		Z6-251 	Z6-271 	Z6-361 	Z6-381 
H:155 mm		Z6-252 	Z6-272 	Z6-362 	Z6-382 
		Z6-253 	Z6-273 	Z6-363 	Z6-383 
		Z6-254 	Z6-274 	Z6-364 	Z6-384 
		Z6-255 	Z6-275 	Z6-365 	Z6-385 
		Z6-256 	Z6-276 	Z6-366 	Z6-386 
		Z6-257 	Z6-277 	Z6-367 	Z6-387 

Ps: All container measurements cover the outer size of the lids.


BOX (perforated)




outside dimension



outside dimension



inside dimension




code	outside dimension	inside dimension	weight
Z6-300	475 x 272 h:100 mm	422 x 258 h:83 mm	1792 g
Z6-301	475 x 272 h:116 mm	422 x 258 h:99 mm	1836 g
Z6-302	475 x 272 h:138 mm	422 x 258 h:121 mm	1980 g
Z6-303	475 x 272 h:185 mm	422 x 258 h:168 mm	2240 g


BOX (non perforated)



outside dimension



outside dimension



inside dimension



code	outside dimension	inside dimension	weight
Z6-305	475 x 272 h:100 mm	422 x 258 h:94 mm	1609 g
Z6-306	475 x 272 h:116 mm	422 x 258 h:110 mm	1714 g
Z6-307	475 x 272 h:138 mm	422 x 258 h:132 mm	1842 g
Z6-308	475 x 272 h:185 mm	422 x 258 h:179 mm	2117 g

PERFORATED LIDS



code	color	weight
Z6-340		1014 g
Z6-341		
Z6-342		
Z6-343		
Z6-344		
Z6-345		

Dimension: 477 x 285 mm

REMOVAL LIDS (non perforated)



code	color	weight
Z6-320		940 g
Z6-321		
Z6-322		
Z6-323		
Z6-324		
Z6-325		

Dimension: 477 x 285 mm

SAFETY LIDS




code	color	weight
Z6-410		903 g
Z6-411		
Z6-412		
Z6-413		
Z6-414		
Z6-415		


Dimension: 483 x 290 mm

BASKET MADE FROM STAINLESS STEEL WIRE MESH

WIRE BASKET



code	outside dimension
Z2-1500	405 x 245 x h:30 mm
Z2-1501	405 x 245 x h:50 mm
Z2-1502	405 x 245 x h:70 mm
Z2-1503	405 x 245 x h:100 mm



WIRE BASKET WITH LID




code	outside dimension
Z3-1030	405 x 245 h:35 mm
Z3-1031	405 x 245 h:55 mm
Z3-1032	405 x 245 h:75 mm
Z3-1033	405 x 245 h:105 mm
Z3-1034	ONLY LID, for baskets




WIRE BASKET WITH PERFORATED METAL PLATE

WIRE BASKET



code	outside dimension
Z2-1501 P	405 x 246 x h:50 mm
Z2-1502 P	405 x 246 x h:70 mm
Z2-1503 P	405 x 246 x h:100 mm



WIRE BASKET WITH LID



code	outside dimension
Z3-1031 P	405 x 246 x h:55 mm
Z3-1032 P	405 x 246 x h:75 mm
Z3-1033 P	405 x 246 x h:105 mm
Z3-1034 P	ONLY LID, for baskets





STERILE CONTAINER SYSTEMS




Standard Model





» *1/2 Half Size*

1/2 STANDARD MODEL

W: 280 mm

L: 310 mm




Specifications		Lid (Perforated)	Lid (Perforated)	Safety Lid	Safety Lid
Height					
		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)	Box (Perforated)
H:110 mm		Z6-500 	Z6-530 	Z6-640 	Z6-670 
		Z6-501 	Z6-531 	Z6-641 	Z6-671 
		Z6-502 	Z6-532 	Z6-642 	Z6-672 
		Z6-503 	Z6-533 	Z6-643 	Z6-673 
		Z6-504 	Z6-534 	Z6-644 	Z6-674 
		Z6-505 	Z6-535 	Z6-645 	Z6-675 
H:130 mm		Z6-506 	Z6-536 	Z6-646 	Z6-676 
		Z6-507 	Z6-537 	Z6-647 	Z6-677 
		Z6-508 	Z6-538 	Z6-648 	Z6-678 
		Z6-509 	Z6-539 	Z6-649 	Z6-679 
		Z6-510 	Z6-540 	Z6-650 	Z6-680 
		Z6-511 	Z6-541 	Z6-651 	Z6-681 
H:150 mm		Z6-512 	Z6-542 	Z6-652 	Z6-682 
		Z6-513 	Z6-543 	Z6-653 	Z6-683 
		Z6-514 	Z6-544 	Z6-654	

BOX (perforated)




outside dimension



outside dimension



inside dimension




code	outside dimension	inside dimension	weight
Z6-600	310 x 272 h:100 mm	253 x 258 h:83 mm	1385 g
Z6-601	310 x 272 h:116 mm	253 x 258 h:99 mm	1417 g
Z6-602	310 x 272 h:138 mm	253 x 258 h:121 mm	1590 g
Z6-603	310 x 272 h:190 mm	253 x 258 h:173 mm	1821 g
Z6-604	310 x 272 h:235 mm	253 x 258 h:218 mm	1947 g


BOX (non perforated)




outside dimension



outside dimension



inside dimension



code	outside dimension	inside dimension	weight
Z6-605	310 x 272 h:100 mm	253 x 258 h:94 mm	1267 g
Z6-606	310 x 272 h:116 mm	253 x 258 h:110 mm	1315 g
Z6-607	310 x 272 h:138 mm	253 x 258 h:132 mm	1427 g
Z6-608	310 x 272 h:190 mm	253 x 258 h:184 mm	1660 g
Z6-609	310 x 272 h:235 mm	253 x 258 h:229 mm	1940 g

PERFORATED LIDS



code	color	weight
Z6-630		727 g
Z6-631		
Z6-632		
Z6-633		
Z6-634		
Z6-635		

Dimension: 312 x 280 mm

REMOVAL LIDS (non perforated)



code	color	weight
Z6-620		629 g
Z6-621		
Z6-622		
Z6-623		
Z6-624		
Z6-625		

Dimension: 312 x 280 mm

SAFETY LIDS



code	color	weight
Z6-700		775 g
Z6-701		
Z6-702		
Z6-703		
Z6-704		
Z6-705		

Dimension: 317 x 285 mm




STERILE CONTAINER SYSTEMS



Standard Model

» *Extra Large Size*

EXTRA LARGE SIZE




Specifications		Lid (Perforated)	
Height			
		Box (Non Perforated)	Box (Perforated)
H:120 mm		Z6-136	Z6-210
		Z6-137	Z6-211
		Z6-138	Z6-212
		Z6-139	Z6-213
		Z6-140	Z6-214
		Z6-141	Z6-215
H:185 mm		Z6-142	Z6-216
		Z6-143	Z6-217
		Z6-144	Z6-218
		Z6-145	Z6-219
		Z6-146	Z6-220
		Z6-147	Z6-221


↻ Ps: All container measurements cover the outer size of the lids.

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
Box and Lid 	Z3-701	539 x 346 h:44 mm
	Z3-702	539 x 346 h:64 mm
	Z3-703	539 x 346 h:84 mm
	Z3-704	539 x 346 h:114 mm
Box 	Z3-705	533 x 340 h:40 mm
	Z3-706	533 x 340 h:60 mm
	Z3-707	533 x 340 h:80 mm
	Z3-708	533 x 340 h:110 mm
Lid 	Z3-709	ONLY LID, for baskets



STAND



code	height
Z3-767	20 mm

STAND



code	height
Z3-768	10 mm

STAND



code	height
Z3-766	5 mm

DIVIDERS



Z2-1731
533 x 337 x 70 mm
Z2-1732
533 x 337 x 120 mm


(-)

DRAPE RETAINER



Z2-1712
533 x 337 x 30 mm

SILICONE MESH MAT



Z2-535
527 x 335 mm

(-)

DIVIDERS



Z2-1721
533 x 337 x 70 mm
Z2-1722
533 x 337 x 120 mm

DRAPE RETAINER



Z2-1711
533 x 337 mm

SILICONE HANDLE



Z2-710
Z2-711
Z2-712
Z2-713
Z2-714
Z2-715

ALUMINIUM BASKET (with lid)



Z3-1119
245 x 130 x 55 mm

ALUMINIUM BASKET (with lid)



Z3-1120
245 x 130 x 70 mm

ALUMINIUM BASKET (with lid)



Z3-1121
245 x 130 x 105 mm

PAPER FILTER




Z2-740
Single use 100 pcs
with indicator. Insidemed
mark.
• Diameter: 184 mm

TEFLON FILTER



Z2-742
Number of sterilization
cycle > 2000
• Diameter: 194 mm

FILTER DISPENSER




Z2-730

IDENTIFICATION LABELS (aluminium)



Z2-700
Z2-701
Z2-702
Z2-703
Z2-704
Z2-705


LABELS (paper)



Z2-673
100 pcs. With indicator.

For small size and
special containers

LABELS (paper)



Z2-672
100 pcs. With indicator.

For full and
Half size containers


(-)

AUTO LOCK MECHANISM



Z1-900
(2 Pcs)

EMULATING INDICATOR CARD (class 6)



Z2-674
Paper 250 Pcs.



SECURITY SEAL (blue)
Z2-663
100 pcs. with indicator.
made of plastic





SECURITY SEAL (blue)
Z2-664
100 pcs. made of plastic






1/1 BIO-BARRIER MODEL

W: 285 mm

L: 600 mm








Specifications ➔	Lid (with barrier)	Lid (with barrier)	Safety Lid Lid (with barrier)	Safety Lid Lid (with barrier)
				
	Box (Non Perforated)	Box (with barrier)	Box (Non Perforated)	Box (with barrier)
	Z7-010 	Z7-040 	Z7-150 	Z7-180 
	Z7-011 	Z7-041 	Z7-151 	Z7-181 
	Z7-012 	Z7-042 	Z7-152 	Z7-182 
	Z7-013 	Z7-043 	Z7-153 	Z7-183 
	Z7-014 	Z7-044 	Z7-154 	Z7-184 
	Z7-015 	Z7-045 	Z7-155 	Z7-185 
	Z7-016 	Z7-046 	Z7-156 	Z7-186 
	Z7-017 	Z7-047 	Z7-157 	Z7-187 
	Z7-018 	Z7-048 	Z7-158 	Z7-188 
	Z7-019 	Z7-049 	Z7-159 	Z7-189 
	Z7-020 	Z7-050 	Z7-160 	Z7-190 
	Z7-021 	Z7-051 	Z7-161 	Z7-191 
	Z7-022 	Z7-052 	Z7-162 	Z7-192 
	Z7-023 	Z7-053 	Z7-	

3/4 BIO-BARRIER MODEL

W: 285 mm

L: 475 mm








Specifications ➔	Lid (with barrier)	Lid (with barrier)	Safety Lid Lid (with barrier)	Safety Lid Lid (with barrier)
				
	Box (Non Perforated)	Box (with barrier)	Box (Non Perforated)	Box (with barrier)
	Z7-240 	Z7-260 	Z7-350 	Z7-370 
	Z7-241 	Z7-261 	Z7-351 	Z7-371 
	Z7-242 	Z7-262 	Z7-352 	Z7-372 
	Z7-243 	Z7-263 	Z7-353 	Z7-373 
	Z7-244 	Z7-264 	Z7-354 	Z7-374 
	Z7-245 	Z7-265 	Z7-355 	Z7-375 
	Z7-246 	Z7-266 	Z7-356 	Z7-376 
	Z7-247 	Z7-267 	Z7-357 	Z7-377 
	Z7-248 	Z7-268 	Z7-358 	Z7-378 
	Z7-249 	Z7-269 	Z7-359 	Z7-379 
	Z7-250 	Z7-270 	Z7-360 	Z7-380 
	Z7-251 	Z7-271 	Z7-361 	Z7-381 
	Z7-252 	Z7-272 	Z7-362 	Z7-382

1/2 BIO-BARRIER MODEL

W: 280 mm

L: 310 mm





Specifications ➔	Lid (with barrier)	Lid (with barrier)	Safety Lid Lid (with barrier)	Safety Lid Lid (with barrier)
				
	Box (Non Perforated)	Box (with barrier)	Box (Non Perforated)	Box (with barrier)
	Z7-500 	Z7-530 	Z7-640 	Z7-670 
	Z7-501 	Z7-531 	Z7-641 	Z7-671 
	Z7-502 	Z7-532 	Z7-642 	Z7-672 
	Z7-503 	Z7-533 	Z7-643 	Z7-673 
	Z7-504 	Z7-534 	Z7-644 	Z7-674 
	Z7-505 	Z7-535 	Z7-645 	Z7-675 
	Z7-506 	Z7-536 	Z7-646 	Z7-676 
	Z7-507 	Z7-537 	Z7-647 	Z7-677 
	Z7-508 	Z7-538 	Z7-648 	Z7-678 
	Z7-509 	Z7-539 	Z7-649 	Z7-679 
	Z7-510 	Z7-540 	Z7-650 	Z7-680 
	Z7-511 	Z7-541 	Z7-651 	Z7-681 
	Z7-512 	Z7-542 	Z7-652 	Z7-682 
	Z7-513			

EXTRA LARGE SIZE



W: 372 mm


L: 600 mm



Specifications ➤	Lid (with barrier)	Lid (with barrier)
	Box (Non Perforated)	Box (with barrier)
	Z7-136	Z7-210
	Z7-137	Z7-211
	Z7-138	Z7-212
	Z7-139	Z7-213
	Z7-140	Z7-214
	Z7-141	Z7-215
	h:130 mm	h:140 mm
	Z7-142	Z7-216
	Z7-143	Z7-217
	Z7-144	Z7-218
	Z7-145	Z7-219
	Z7-146	Z7-220
	Z7-147	Z7-221
	h:195 mm	h:205 mm

➤ Ps: All container measurements cover the outer size of the lids.

SILICONE MAT


**Z2-526**

520 x 230 mm

Z2-528

475 x 230 mm

SILICONE MESH MAT

**Z2-531**

485 x 250 mm

SILICONE HANDLE

**Z2-710****Z2-711****Z2-712****Z2-713****Z2-714****Z2-715**

PAPER FILTER

**Z2-740**


Single use 100 pcs
with indicator. Insidemed
mark. Diameter: 184 mm

ALUMINIUM FILTER RETAINER


**Z2-722**

Diameter: 189 mm


FILTER DISPENSER

**Z2-730**

IDENTIFICATION LABELS (aluminium)

**Z2-700****Z2-701****Z2-702****Z2-703****Z2-704****Z2-705**


LABELS (paper)

**Z2-673**

100 pcs. With indicator.

For small size and
special containers

LABELS (paper)

**Z2-672**

100 pcs. With indicator.

For full and
Half size containers

SECURITY SEAL (blue)

**Z2-663**

100 pcs. with indicator.
made of plastic

SECURITY SEAL (blue)

**Z2-664**

100 pcs. made of plastic

SECURITY SEAL (blue)

**Z2-661**

100 pcs. with number.
made of plastic

SECURITY SEAL (red)

**Z2-662**

100 pcs. with indicator.
made of plastic

SECURITY SEAL (red)

**Z2-660**


100 pcs. made of plastic

SECURITY SEAL (red)

**Z2-665**


100 pcs. with number.
made of plastic

SILICONE MAT



Z2-524
380 x 230 mm
Z2-525
440 x 230 mm

SILICONE MESH MAT



Z2-531
485 x 250 mm

SILICONE HANDLE



Z2-710	Red
Z2-711	Yellow
Z2-712	Green
Z2-713	Blue
Z2-714	Grey
Z2-715	Black

PAPER FILTER




Z2-740
Single use 100 pcs
with indicator. Insidemed
mark.
• Diameter: 184 mm

ALUMINIUM FILTER RETAINER




Z2-722
Diameter: 189 mm

FILTER DISPENSER




Z2-730

IDENTIFICATION LABELS (aluminium)



Z2-700	Red
Z2-701	Blue
Z2-702	Green
Z2-703	Yellow
Z2-704	Grey
Z2-705	Black


LABELS (paper)



Z2-673
100 pcs. With indicator.

For small size and
special containers

LABELS (paper)



Z2-672
100 pcs. With indicator.

For full and
Half size containers

SECURITY SEAL (blue)



Z2-663
100 pcs. with indicator.
made of plastic

SECURITY SEAL (blue)



Z2-664
100 pcs. made of plastic

SECURITY SEAL (blue)



Z2-661
100 pcs. with number.
made of plastic

SECURITY SEAL (red)



Z2-662
100 pcs. with indicator.
made of plastic

SECURITY SEAL (red)



Z2-660
100 pcs. made of plastic

SECURITY SEAL (red)




Z2-665
100 pcs. with number.
made of plastic

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
	Z3-752	244 x 253 h:40 mm
	Z3-753	244 x 253 h:60 mm
	Z3-754	244 x 253 h:80 mm
	Z3-755	244 x 253 h:110 mm

STAND



code	height
Z3-767	20 mm




STAND



code	height
Z3-768	10 mm



STAND



code	height
Z3-766	5 mm



BASKET MADE FROM STAINLESS STEEL WIRE MESH



code	outside dimension
Z2-490	246 x 246 x h:30 mm
Z2-491	246 x 246 x h:50 mm
Z2-492	246 x 246 x h:70 mm
Z2-493	246 x 246 x h:100 mm

