

VATS/MICS



NERVE AND VESSEL HOOKS

	ITEM NO.	LENGTH (MM)	DESCRIPTION
	GV105.30	300	NERVE AND VESSEL HOOK, ANGLED 90°, LENGTH 6 MM
	GV104.30	300	NERVE AND VESSEL HOOK, ANGLED 90°, LENGTH, 6 MM, BALL POINT
	GV106.30	300	NERVE AND VESSEL HOOK, ANGLED 90°, LENGTH 10 MM



KNOT PUSHER

	ITEM NO.	LENGTH (MM)	Ø	DESCRIPTION
	GV065.25	250	5MM	STRAIGHT



MAGNET RETRIEVER

	ITEM NO.	LENGTH (MM)	Ø	DESCRIPTION
	GV071.25	250	4MM	STRAIGHT

VATS/MICS

SUTURE CATCHER



	ITEM NO.	LENGTH (MM)	Ø	DESCRIPTION
	GV070.30	230	2.0	STRAIGHT



DISSECTION INSTRUMENT CLAMP

ITEM NO.	LENGTH (MM)	Ø	DESCRIPTION
GV575.25	25MM	5.0	SPREADING INSTRUMENT , JAW FINE SERRATED




MEASURING INSTRUMENT

ITEM NO.	LENGTH (MM)	Ø	DESCRIPTION
GV035.28	280MM	3.5	MEASURING INSTRUMENT, SIZES, 0-60 MM

FORCEPS

D'AMICO BIOPSY FORCEPS



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV302.13	6 MM WIDTH OVAL RING JAWS CURVED, WITH RATCHET	5MM TAPERED TO 7MM	9" (23CM)	13" (33CM)
GV303.13	6 MM WIDTH OVAL RING JAWS CURVED WITHOUT RATCHET			



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV312.13	7 MM WIDTH OVAL CUPPED JAWS CURVED WITH RATCHET	5MM TAPERED TO 7MM	9" (23CM)	13" (33CM)
GV313.13	7 MM WIDTH OVAL CUPPED JAWS CURVED WITHOUT RATCHET			

FORCEPS

D'AMICO MEDIASTINOSCOPY FORCEPS



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV300.15	6 MM WIDTH OVAL RING JAWS CURVED, WITH RATCHET	7MM	11" (28CM)	15" (38CM)
GV301.15	6 MM WIDTH OVAL RING JAWS CURVED WITHOUT RATCHET			



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV310.15	7 MM WIDTH OVAL CUPPED JAWS CURVED WITH RATCHET	7MM	11" (28CM)	15" (38CM)
GV311.15	7 MM WIDTH OVAL CUPPED JAWS CURVED WITHOUT RATCHET			

NODULE CLAMPS

ROCCO NODULE CLAMPS



GV946.34




GV947.34

ITEM NO.	DESCRIPTION	SUGGESTED NODULE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV946.34	CURVED LEFT RING JAWS, 3CM INNER DIAMETER AND 3.7 CM OUTER DIAMETER WITH RATCHET	UP TO 1CM	10MM	9½" (24CM)	13¼" (33.5CM)
GV947.34	CURVED LEFT RING JAWS 4.5 CM INNER DIAMETER AND 5.2 CM OUTER DIAMETER, WITH RATCHET	BETWEEN 1 & 1.5CM			


GRASPING CLAMPS

GONZALEZ-RIVAS LUNG GRASPING CLAMP



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV230.25	HIGH TENSION CURVED LEFT 20 MM OVAL RING JAWS WITH RATCHET	7MM TAPERED TO 10MM	5¾" (14.5CM)	9¾" (24.5 CM)

KHAN DECORTICATION GRASPER



KAHN DECORTICATION GRASPER JAWS COMPARED TO
FOERSTER CLAMP JAWS

GV940.00—KHAN DECORTIFICATION
GRASPER

GV920.00—FOERSTER LUNG
CLAMP



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV940.34	ENHANCED SERRATIONS CURVED LEFT JAWS WITH RATCHET	10MM	9½" (24 CM)	13¼" (33.5 CM)

GRASPING CLAMPS

FOERSTER LUNG GRASPING CLAMP



#GV229.34



#GV237.42

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV229.34	CURVED LEFT 20 MM OVAL RING JAWS WITH RATCHET	5MM TAPERED TO 7MM	9¼" (23.5CM)	13¼" (33.5CM)
GV237.42			12½" (32CM)	16½" (41.9 CM)

GRASPING CLAMPS

FOERSTER LUNG GRASPING CLAMP




#GV227.29



#GV228.29

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV227.29	CURVED RIGHT JAWS (LEFT-HAND STYLE) 12 MM OVAL RING JAWS, WITH RATCHET	7MM	7½" (19CM)	11½" (29CM)
GV228.29	CURVED LEFT JAWS 12 MM OVAL RING JAWS WITH RATCHET			




#GV226.34

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV226.34	CURVED LEFT JAWS 12 MM OVAL RING JAWS WITH RATCHET	7 MM	9¼" (23.5 CM)	13¼" (33.5 CM)

GRASPING CLAMPS


FOERSTER LUNG GRASPING CLAMP



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV921.29	CURVED RIGHT JAWS (LEFT-HAND STYLE) 20 MM OVAL RING JAWS WITH RATCHET	10 MM	7½" (19CM)	11½" (29CM)
GV922.29	CURVED LEFT JAWS 20 MM OVAL RING JAWS WITH RATCHET			

GRASPING CLAMPS

FOERSTER LUNG GRASPING CLAMP



#GV923.33



#GV919.43

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV923.33	DEEP CURVED LEFT JAWS 20 MM OVAL RING JAWS WITH RATCHET	10 MM	9½" (24 CM)	11½" (29CM)
GV919.43	CURVED LEFT JAWS 20 MM OVAL RING JAWS WITH RATCHET		12½" (32CM)	16¾" (42.4 CM)

GRASPING CLAMPS


FOERSTER LUNG GRASPING CLAMP



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV920.34	CURVED LEFT JAWS 20 MM OVAL RING JAWS	10 MM	9½" (24 CM)	13¼" (33.5 CM)
GV925.43	WITH RATCHET		12½" (32 CM)	16¾" (42.4 CM)
GV924.34	CURVED LEFT JAWS 25 X 35 MM OVAL RING JAWS WITH RATCHET		9½" (24 CM)	13¼" (33.5 CM)

GRASPING CLAMPS


DUVAL LUNG GRASPING CLAMPS



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV926.34	CURVED LEFT JAWS	10 MM	9½" (23.5 CM)	13½" (34 CM)
GV928.43	1.0 CM WIDTH TRIANGULAR RING JAWS WITH RATCHET		12½" (32 CM)	16¾" (42.4 CM)
GV927.34	CURVED LEFT JAWS		9¼" (23.5 CM)	13½" (34 CM)
	2.5 CM WIDTH TRIANGULAR RING JAWS WITH RATCHET			

GRASPING CLAMPS

NODE GRASPING CLAMPS



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV222.33	STRAIGHT JAWS 11 MM OBLONG RING JAWS WITH RATCHET	5MM	9" (23 CM)	13" (33CM)
GV223.33	CURVED JAWS 11 MM OBLONG RING JAWS WITH RATCHET			
GV225.13	CURVED JAWS 11 MM OBLONG RING JAWS WITHOUT RATCHET			
GV220.43	CURVED JAWS 11 MM OBLONG RING JAWS WITH RATCHET	5 MM TAPERED TO 7 MM	12½" (32 CM)	16¾" (42.4 CM)
GV221.43	CURVED JAWS 11 MM OBLONG RING JAWS WITHOUT RATCHET			

DISSECTORS

GONZALEZ-RIVAS DISSECTOR



GV.210
GV.211
GV.214
GV.215
GV.216
GV.217



GV.212
GV.213

GONZALEZ-RIVAS DISSECTOR REGULAR JAW TIPS
COMPARED TO FINE JAW TIPS



GV210.19 &
GV211.29

REGULAR JAW TIPS




GV212.29 &
GV213.29

REGULAR JAW TIPS

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV210.29	MIXTER-STYLE 1 X 2 DeBAKEY JAWS WITH RATCHET	5 MM	7½" (19 CM)	11½" (29 CM)
GV211.29	MIXTER-STYLE 1 X 2 DeBAKEY JAWS WITHOUT RATCHET			
GV212.29	MIXTER-STYLE 1 X 2 DeBAKEY JAWS FINE TIPS WITH RATCHET			
GV213.29	MIXTER-STYLE 1 X 2 DeBAKEY JAWS FINE TIPS WITHOUT RATCHET			
GV214.34	MIXTER-STYLE 1 X 2 DeBAKEY JAWS WITH RATCHET	5 MM TO 7 MM TAPERED	9½" (24 CM)	13½" (34 CM)
GV215.34	MIXTER-STYLE 1 X 2 DeBAKEY JAWS WITHOUT RATCHET			
GV216.42	MIXTER-STYLE 1 X 2 DeBAKEY JAWS WITH RATCHET	5 MM TO 7 MM TAPERED	12½" (32 CM)	16½" (41.9 CM)
GV217.42	MIXTER-STYLE 1 X 2 DeBAKEY JAWS WITHOUT RATCHET			


DISSECTORS

SNAKE DISSECTOR



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV284.34	CURVED DISSECTOR JAWS WITH RATCHET	5 MM TAPERED TO 7 MM	9¼" (23.5 CM)	13½" (34.4 CM)
GV285.34	CURVED DISSECTOR JAWS WITHOUT RATCHET			
GV286.43	CURVED DISSECTOR JAWS WITH RATCHET		12½" (32 CM)	16¾" (42.4 CM)
GV287.43	CURVED DISSECTOR JAWS WITHOUT RATCHET			


DENNIS DISSECTOR



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV280.34	CURVED LEFT JAWS SERRATED JAWS WITH RATCHET	7 MM	9" (23 CM)	13¼" (33.5 CM)
GV281.34	CURVED LEFT JAWS SERRATED JAWS WITHOUT RATCHET			
GV282.43	CURVED LEFT JAWS SERRATED JAWS WITH RATCHET	5 MM TAPERED TO 7 MM	12.5" (32 CM)	16¾" (42.4 CM)
GV283.43	CURVED LEFT JAWS SERRATED JAWS WITHOUT RATCHET			

CLAMPS

DEBAKEY COOLEY CLAMPS



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV202.34	DEBAKEY COOLEY JAWS WITH RATCHET	10 MM	9½" (24 CM)	13½" (34 CM)
GV203.34	DEBAKEY COOLEY JAWS WITHOUT RATCHET			

CLAMPS


HARKEN CLAMPS



ITEM NO.	DESCRIPTION		SHAFT	WORKING LENGTH	TOTAL LENGTH
GV250.29	HARKEN #1	7.5 CM JAWS TAPERED 1x2 DEBAKEY JAWS WITHOUT RATCHET	5 MM TO 7 MM TAPERED	7½" (19 CM)	11½" (29 CM)
GV251.29		7.5 CM JAWS TAPERED 1x2 DEBAKEY JAWS WITH RATCHET			
GV252.29	HARKEN #2	7.5 CM JAWS TAPERED 1x2 DEBAKEY JAWS WITHOUT RATCHET			
GV253.29		7.5 CM JAWS TAPERED 1x2 DEBAKEY JAWS WITH RATCHET			
GV224.36	HARKEN #1	9.5 CM JAWS TAPERED 1x2 DEBAKEY JAWS WITHOUT RATCHET	10 MM	9⅞" (25 CM)	14" (35.6 CM)
GV225.36		9.5 CM JAWS TAPERED 1x2 DEBAKEY JAWS WITH RATCHET			
GV226.36	HARKEN #2	9.5 CM JAWS TAPERED 1x2 DEBAKEY JAWS WITHOUT RATCHET			
GV227.36		9.5 CM JAWS TAPERED 1x2 DEBAKEY JAWS WITH RATCHET			

CLAMPS


ALLIS FORCEPS (CLAMP)



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV270.34	CURVED LEFT JAWS 6 MM WIDTH JAWS WITH RATCHET	7 MM	9" (23 CM)	13¼" (33.5CM)
GV274.29	STRAIGHT JAWS 4 MM WIDTH JAWS WITH RATCHET	5 MM TAPERED TO 7 MM	7½" (19 CM)	11½" (29 CM)
GV275.29	CURVED LEFT JAWS 4 MM WIDTH JAWS WITH RATCHET			

CLAMPS

D'AMICO DeBAKEY CLAMPS



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV219.37	CURVED LEFT JAWS TAPERED 1x2 DeBAKEY JAWS 11CM JAW LENGTH WITHOUT RATCHET	10 MM	10" (25.4)	14½" (37 CM)
GV221.37	CURVED LEFT JAWS TAPERED 1x2 DeBAKEY JAWS 11CM JAW LENGTH WITH RATCHET			
GV222.38	STRAIGHT JAWS TAPERED 1x2 DeBAKEY JAWS 11.5 CM JAW LENGTH WITH RATCHET		10½" (27 CM)	14¾" (37.5 CM)

CHITWOOD DeBAKEY CLAMPS

GV240.31
GV241.31



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV240.31	6 CM JAW LENGTH CURVED LEFT JAWS 1 x 2 DeBAKEY JAWS WITH RATCHET	5 MM	8¼" (21 CM)	12¼" (31 CM)
GV241.31	6 CM JAW LENGTH CURVED LEFT JAWS 1 x 2 DeBAKEY JAWS WITHOUT RATCHET			

CLAMPS

CHITWOOD DeBAKEY CLAMPS



GV912.33



GV912.36



GV232.36



GV232.36

2 x 3 DeBAKEY JAWS




GV912.36

1 x 2 DeBAKEY JAWS

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV912.33	6.5 CM JAW LENGTH 1x2 DeBAKEY JAWS CURVED LEFT JAWS WITH RATCHET	10 MM	9" (23 CM)	12¾" (32.5 CM)
GV912.35	10 CM JAW LENGTH 1x2 DeBAKEY JAWS CURVED LEFT JAWS WITH RATCHET		10" (25.4 CM)	14¼" (36 CM)
GV232.36	10 CM JAW LENGTH 2x3 DeBAKEY JAWS CURVED LEFT JAWS WITH RATCHET			

CLAMPS

CHITWOOD DEBAKEY CLAMPS



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV914.27	7 CM JAW LENGTH CURVED LEFT JAW 1x2 DEBAKEY JAWS WITH RATCHET	10 MM	6 ³ / ₄ " (17 CM)	10 3/8" (26.4 CM)
GV915.27	7 CM JAW LENGTH CURVED RIGHT JAW 1x2 DEBAKEY JAWS WITH RATCHET			

CLAMPS

CHITWOOD DeBAKEY CLAMPS




ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV234.21	8 CM JAW LENGTH CURVED LEFT JAWS 2 X 3 DeBAKEY JAWS ANGLED SHANKS	10 MM	5 $\frac{1}{8}$ " (13CM)	8" (20.5 CM)



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV913.25	7 CM JAW LENGTH CURVED LEFT JAWS 1 X 2 DeBAKEY JAWS ANGLED SHANKS	10 MM	6 $\frac{3}{4}$ " (17 CM)	9 $\frac{3}{4}$ " (24.8 CM)

SCISSORS

SUPER CUT SCISSORS



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV716.27	STRAIGHT FINE BLADES 1 CM BLADE LENGTH BLUNT	5 MM	6½" (16.5 CM)	10½" (27 CM)
GV718.27	CURVED FINE BLADES 1 CM BLADE LENGTH BLUNT		6½" (16.5 CM)	10½" (27 CM)
GV738.29	CHITWOOD DISSECTING SCISSORS CURVED BLADES 3 CM BLADE LENGTH BLUNT TIPS		7¼" (18.5 CM)	11½" (29 CM)

SCISSORS

SUPER CUT SCISSORS



METZENBAUM SCISSORS

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV710.42	CURVED BLADES 3.1 CM BLADE LENGTH BLUNT TIPS	10 MM	12½” (32 CM	16½” (42.1CM)
GV910.29	CURVED BLADES 3.5 CM BLADE LENGTH BLUNT TIPS		7¾” (19.5 CM)	11½” (29 CM)
GV911.29	CURVED BLADES 3.5 CM BLADE LENGTH BLUNT TIPS			



GV910.29



GV911.29



SURGICAL VATS/MIS

ROCCO ADVENTITIA ELEVATOR (FORCEPS)



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV500.33	DIAMOND DUST JAWS CURVED DOWN WITHOUT RATCHET RING HANDLE CURVED DOWN	5 MM	8" (20 CM)	12¾" (32.5 CM)

NEEDLE HOLDERS

SPRING STYLE NEEDLE HOLDERS



ITEM NO.	DESCRIPTION	SUTURE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV648.39	STRAIGHT DIAMOND DUST JAWS WITH RATCHET	5-0 AND SMALLER	5MM	6½" (16.5CM)	153/8" (38.5CM)
GV650.39	CURVED DIAMOND DUST JAWS WITH RATCHET				

RING HANDLE NEEDLE HOLDERS




ITEM NO.	DESCRIPTION	SUTURE	SHAFT	WORKING	TOTAL LENGTH
GV618.27	STRAIGHT JAWS FINE DIAMOND DUST JAWS WITH RATCHET	5-0 AND SMALLER	5 MM	6½" (16.5CM)	10½" (27CM)

NEEDLE HOLDERS

RING HANDLE NEEDLE HOLDERS



ITEM NO.	DESCRIPTION	SUTURE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV636.29	STRAIGHT JAWS CARBON-INLAY JAWS GOLD RING HANDLE WITH RATCHET	4-0 AND SMALLER	7MM	6¼" (16CM)	11¼" (28.5CM)
GV637.29	CURVED JAWS CARBON-INLAY JAWS GOLD RING				
GV638.29	STRAIGHT JAWS CARBON-INLAY JAWS GOLD RING	2-0 AND SMALLER	7MM	6¼" (16CM)	11¼" (28.5CM)
GV639.29	CURVED JAWS CARBON-INLAY JAWS GOLD RING				




ITEM NO.	DESCRIPTION	SUTURE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV644.45	STRAIGHT JAWS CARBON-INLAY JAWS GOLD RING HANDLE WITH RATCHET	2-0 AND SMALLER	7 MM	12½" (32CM)	17½" (44.5CM)
GV645.45	CURVED JAWS CARBON-INLAY JAWS GOLD RING HANDLE WITH RATCHET				

NEEDLE HOLDERS

RING HANDLE NEEDLE HOLDERS



ITEM NO.	DESCRIPTION	SUTURE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV610.34	STRAIGHT JAWS CARBON-INLAY JAWS GOLD RING HANDLE WITH RATCHET	2-0 AND SMALLER	10MM	8¼" (21CM)	13½" (34CM)
GV612.34	CURVED JAWS CARBON-INLAY JAWS GOLD RING HANDLE WITH RATCHET				



ITEM NO.	DESCRIPTION	SUTURE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV640.34	STRAIGHT JAWS CARBON-INLAY JAWS GOLD RING HANDLE WITH RATCHET	2-0 AND SMALLER	7MM	8¼" (21CM)	13½" (34CM)
GV642.34	CURVED JAWS CARBON-INLAY JAWS GOLD RING HANDLE WITH RATCHET				

KNOT PUSHERS

GONZALEZ-RIVAS KNOT PUSHER



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV826.38	CURVED UP	5MM	11½" (27CM)	14¾" (37.5CM)
GV826.43			12½" (32CM)	16¾" (42.4CM)

CHITWOOD KNOT PUSHER



ITEM NO.	DESCRIPTION	SUTURE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV827.28	"TRAP MODEL"	2-0 & SMALLER	2-0 & SMALLER	5⅛" (13CM)	11" (28CM)
GV827.31				5½" (14CM)	12" (30.5CM)
GV827.35				6¾" (17CM)	13⅝" (34.6CM)
GV827.39				8½" (21.5CM)	15 3/8" (39.1CM)


KNOT PUSHERS

CHITWOOD KNOT PUSHER



ITEM NO.	DESCRIPTION	SUTURE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV282.32	ORIGINAL "TRAP MODEL"	2-0 & SMALLER	5MM	6 $\frac{1}{8}$ " (15.5CM)	12 $\frac{1}{2}$ " (32CM)
GV282.39	ORIGINAL "TRAP MODEL"			8 $\frac{1}{2}$ " (21.5CM)	15 $\frac{3}{8}$ " (39.1CM)
GV828.40	ORIGINAL "TRAP MODEL" LONG HANDLE			7 $\frac{1}{4}$ " (18.5CM)	15 $\frac{3}{4}$ " (40CM)

CHITWOOD SUPER CUT SUTURE CUTTER



ITEM NO.	DESCRIPTION	SUTURE SIZE	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV838.28	"TRAP MODEL"	2-0 & SMALLER	5MM	7" (18CM)	11" (28CM)
GV840.33				9" (23CM)	13" (33CM)
GV842.38				11" (28CM)	15" (38CM)

SUCTION INSTRUMENTS

WOLF



VARIABLE SUCTION REGULATOR

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV908.31	CURVED BLUNT TIP 6 MM WIDTH TIP	3.2MM	6¾" (17CM)	12¼" (31CM)
GV909.31	VARIABLE SUCTION CURVED BLUNT TIP 6 MM WIDTH TIP			
GV910.36	CURVED BLUNT TIP 6 MM WIDTH TIP		8¾" (22CM)	14¼" (36CM)
GV911.36	VARIABLE SUCTION CURVED BLUNT TIP 6 MM WIDTH TIP			
GV912.36	CURVED BLUNT TIP 6 MM WIDTH TIP	32.5 TO 5MM	8¾" (22CM)	14¼" (36CM)
GV913.36	VARIABLE SUCTION CURVED BLUNT TIP 6 MM WIDTH TIP			

SUCTION INSTRUMENTS

DENNIS



VARIABLE SUCTION REGULATOR

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV914.36	CURVED BLUNT TIP 7.5 MM WIDTH TIP	5MM	8 ³ / ₄ " (22CM)	14 ¹ / ₄ " (36CM)
GV915.36	VARIABLE SUCTION CURVED BLUNT TIP 7.5 MM WIDTH TIP			
GV916.46	CURVED BLUNT TIP 7.5 MM WIDTH TIP		12 ¹ / ₂ " (32CM)	17 ³ / ₄ " (45.5CM)
GV917.46	VARIABLE SUCTION CURVED BLUNT TIP 7.5 MM WIDTH TIP			

SUCTION INSTRUMENTS

D'AMICO



VARIABLE SUCTION REGULATOR



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV900.46	STRAIGHT BLUNT TIP 5 MM WIDTH TIP	5MM	12½" (32CM)	17¾" (45.5CM)
GV901.46	VARIABLE SUCTION STRAIGHT BLUNT TIP 5 MM WIDTH TIP			

SUCTION ACCESSORIES



ITEM NO.	DESCRIPTION	TOTAL LENGTH
GV905.49	BRUSH SET, PACK OF 10, FOR 5 MM D'AMICO AND 3.2 MM WOLF AND DENNIS SUCTION INSTRUMENTS	19" (48.3CM)
GV922.49	BRUSH SET, PACK OF 10, FOR 10 MM D'AMICO SUCTION INSTRUMENTS	

SUCTION INSTRUMENTS

D'AMICO



VARIABLE SUCTION REGULATOR

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV920.46	STRAIGHT BLUNT TIP 10 MM WIDTH TIP	10MM	12½" (32CM)	17¾" (45.5CM)
GV921.46	VARIABLE SUCTION STRAIGHT BLUNT TIP 10 MM WIDTH TIP			
GV924.46	CURVED BLUNT TIP 10 MM WIDTH TIP			
GV925.46	VARIABLE SUCTION CURVED BLUNT TIP 10 MM WIDTH TIP			

DENNIS RIB CUTTER



ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV900.34	10 MM BLADE 16 MM CUTTING WIDTH STRAIGHT SHAFT	10MM	11½" (29CM)	13½" (34CM)
GV902.34	15 MM BLADE 25 MM CUTTING WIDTH STRAIGHT SHAFT			
GV904.34	20 MM BLADE 25 MM CUTTING WIDTH STRAIGHT SHAFT			




GV900.34



GV902.34



GV904.34



ITEM NO.	BLADE HEIGHT	CUTTING WIDTH
GV900.34	10MM	16MM
GV902.34	15MM	25MM
GV904.34	20MM	25MM

CHITWOOD SUTURE PULLER



GV816.29



GV820.29

ITEM NO.	DESCRIPTION	SHAFT	WORKING LENGTH	TOTAL LENGTH
GV816.29	STRAIGHT 1.75MM TIP	3.1MM	7½" (19CM)	11½" (29CM)
GV820.29	STRAIGHT 3.1MM TIP			